




System Integration


Part 1. Integrated system design

Part 2. Integrating individual components of the system



Overall hardware block diagram I



Overall hardware block diagram II





Overall hardware block diagram III



Part 1. Integrated system design

- Framework flow
- Iterative flow
- Back to user
- Parameters flow





Neuroimaging system

Imaging field of view: cylinder length $200 \times 200 \times 200 \text{ mm}^3$

Need a spatial resolution of $3 \times 3 \times 3 \text{ mm}$

Estimate that we want a magnetic field strength of at least 45 mT

Maximum length of the magnet is 500 mm to position the brain at the centre

So the RF coil is specified as cylinder $240 \times 240 \times 240 \text{ mm}$


Gradient coil length set at 500 mm

Diameter 300 mm (to avoid reducing Tx/Rx efficiency)


Effect of RF shield on RF coil performance

A

Dome-helix




Elliptical solenoid




B


Domehelix
Cylindrical shield




Elliptical solenoid
Cylindrical shield




Dome-helix
Elliptical shield




Cylindrical shield




Elliptical shield



Decrease in SNR as function of shield-to-coil distance






Optimization via many different methods

Most methods model individual magnets as dipoles

Start with a reasonable guess in terms of diameter, length, number of rings based on previous designs (although not critical)


In this case, 25 rings, 20 mm in-between rings



Can use different approaches, target field, boundary element methods etc.....

Gradient design algorithm:

1. Initialization:
 - a. Obtain gradient coil length and radius from magnet dimensions
 - b. Specify target gradient field vector \mathbf{k}
 - c. Specify the wire thickness, and minimum wire spacing
 - d. Specify maximum linearity
2. Find expansion coefficients \sim for a range of the regularization parameter λ
3. Determine gradient efficiencies and linearities for all values of λ
4. Find most efficiency coil corresponding to the maximum allowable linearity error.



Consider the effect of B_0 inhomogeneity and gradient strength (bandwidth)


$$\frac{1}{T_2^*} = \frac{1}{T_2} + \gamma \Delta B_0$$

$$\gamma \Delta B_0 \gg \frac{1}{T_2} \quad T_2^* \approx \frac{1}{\gamma \Delta B_0}$$

$BW = FOV \times \text{gradient strength}$

$$\Delta z_{eff} = N_z \Delta z \frac{\tanh\left(\frac{\pi FWHM}{2BW}\right)}{1 - e^{-\frac{\pi FWHM N_z}{2BW}}}$$

where BW is the maximum obtainable readout bandwidth of the gradient system determined with the least efficient gradient, Δz the digital resolution and N_z the number of points.



Iterative loop

FOV 200 x 200 x 200, spatial resolution 3 x 3 x 3

Data matrix: 80 x 66 x 66 assuming 0.5 mm frequency encoding blurring

Magnet: 25 rings, 20 mm ring spacing

Optimization of spacing and diameter gives 4329 ppm

Gradient coils: optimization gives an X-gradient with efficiency 0.24 mT/m/A. Maximum 30 Amps gives an effective spatial resolution of ~10 mm.



Magnet: add one ring to increase B_0 homogeneity

Optimization of spacing and diameter gives 1075 ppm

Gradient length increases, 0.26 mT/m/A, resolution ~4.1 mm



Magnet: add one ring to increase B_0 homogeneity


Optimization of spacing and diameter gives 381 ppm

Gradient length increases, resolution ~2.6 mm

Magnet: increase target B_0 in 1 mT steps
Change length and diameter, while
Maintaining target resolution



Check that the RF side is adequate



Solenoid: 240 x 240 x 240 mm

B_1 uniformity > 20%

Transmit efficiency 26 microTesla per $\sqrt{\text{Watt}}$ input power

1 mW maximum output power of spectrometer

RF amplifier gain 54 dB

Maximum flip angle 180°

Pulse duration 50 microseconds

Maximum output power is 9 Watts

RF pulse excites bandwidth of ~ 20 kHz which is greater than ΔB_0

Part 2. Integrating individual components of the system

RF amplifier


What is maximum input to RF amplifier?

What is the maximum duty cycle/pulse length the amplifier can deal with?

What is the shortest pulse that it can amplify?

What is the dynamic range of the amplifier (soft pulses)?

How much reflected power can it take?



Assume that we have an RF coil that is impedance matched at ~50 Ohms

Characterizing the RF amplifier performance




TABLE 3 | Measured pulse characteristics at 2.15 MHz at the output of the RF amplifier.


Blanked noise voltage	0.0066 ± 0.0021 V
Bias enable transmit voltage	1.84 V
Unblanked noise voltage	0.0096 ± 0.004 V
Pulse rise transition duration (10–90% of peak output voltage)	<1 μ s
Rising pulse overshoot	<1%
Pulse overshoot ringing/decay time	<1 μ s
Pulse tilt positive/negative	<1%
Pulse falling transition duration (90–10% of peak output voltage)	<1 μ s

Characterizing the RF amplifier performance

Output of
RF amplifier



After Tx/Rx
switch




Pick-up coil loosely
coupled to RF coil




Characterizing the RF amplifier performance

voltage gain (dB)



voltage gain (dB) @ 2.15 MHz



Looking inside the RF amplifier







Figure 6.17 Schematic of a broadband matching network.

Design of a simple ~10 Watt output amplifier



Gradient amplifier



Simple gradient coil characteristics

Gradient coil has a resistance, inductance (L), capacitance C from the different turns, frequency dependent inductive reactance (ωL), and frequency dependent capacitive reactance ($1/\omega C$)


So there is a resonant frequency (or several resonant frequencies) from the gradient coils which we have to be aware of since they typically occur in the \sim MHz range, good shielding between gradients and RF is required. Otherwise they will couple together VERY strongly.

The impedance of the gradient coils is frequency dependent, so their characteristics depend on whether the gradient is being turned on (high frequency components) or is in steady-state (very low frequency)

Typical numbers: inductance 200 microhenries, resistance 0.4 Ohms.

At 1 kHz, $\omega L = 1.2$ Ohms, at 10 kHz $\omega L = 12$ Ohms. So complex impedance is $0.4 + j1.2$ or $0.4 + j12$

Feedforward and feedback controller for “correcting” the gradient waveform



When running a constant-current control loop that drives a voltage across an inductive load, the load's impedance directly affects the control loop response and stability, and so to be able to make the amplifier work well with a variety of gradient coils (a range of both load inductances and resistances), the control loop values need to be adjusted.

Technical Data Draft

Gradient Amplifier GA-30-15

Part-No. 410000

Output Channels 1

Output Current ± 30 A peak / ± 20 A average

Output Voltage ± 15 V

Load Impedance 0.1 – 4.8 Ohms
0 – 20 mHenry

Digital Resolution 16 bit internal DAC / 1 mA steps

Input Voltage for channel 0 ± 10 V



System specifications

Supply voltage 100 – 400 Vdc
Input current max 30 A

Load inductance 30 μ H – 5 mH
Load resistance 0 Ω - 1 Ω


Environmental requirements:

Ambient temperature 10 °C to 40 °C
Ambient humidity 30 to 75 % non-condensing
Transport temperature -25 °C to +60 °C
Cooling Air cooling (front in, rear out)



3-Axis MRI
Gradient
Amplifier
175A 350V


Rise time of the gradient coils for a step input function is given by L/R




Match the DC resistance of gradient coil to designed load resistance of the amplifier achieves maximum continuous gradient.

Reducing the impedance of the gradient coil decreases the rise time, but reduces peak gradient

underdamped
response



corrected



Illustrative examples of passive “impedance matching” for a Bruker amplifier

Table 7.1. Micro5 Probe and Gradient Specifications

Gradients	XYZ
Gradient strength	4.8 G/cm/A
ID/OD	19/40 mm
Linearity +-1.3% peak-peak +- 1.6% peak-peak +- 2.1% peak-peak	18 mm sphere 19 mm sphere 20 mm sphere
Inductance	10 - 20 μ H
Resistance	=< 120 m Ω
Rise time, 0-40A, 120V	< 50 μ s
Cooling	air or water
Maximum current tested	40 A

Table 7.2. B-AFPA-40 Dip Switch Setting

Gradient system	Gradient	C Dip Switch	R Dip Switch	C	R
Micro5	X	0000 1000	0001 0000	1.5 nF	20 k Ω
Micro5	Y	0000 0100	0001 0000	1.0 nF	20 k Ω
Micro5	Z	0000 1000	0000 1000	1.5 nF	15 k Ω

Table 7.4. Micro2.5 probe and gradient specifications

Gradients	XYZ
Gradient strength	2.5 G/cm/A
ID/OD	40/72 mm
Linearity +-1.8% peak-peak +- 2.2% peak-peak +- 3.0% peak-peak	36 mm sphere 38 mm sphere 40 mm sphere
Inductance	=< 100 μ H
Resistance	=< 400 m Ω
Rise time, 0-40A, 120V	< 110 μ s
Cooling	air or water
Maximum current tested	40 A

Gradient system	Gradient	C Dip Switch	R Dip Switch	C	R
Micro2.5	X	00001100	00010000	2.5 nF	20 k Ω
Micro2.5	Y	00001100	00010000	2.5 nF	20 k Ω
Micro2.5	Z	00001100	00010000	2.5 nF	20 k Ω



Table 7.8. Mini05 Probe and Gradient Specifications

gradients	XYZ
gradient strength	0.5 G/cm/A
ID/OD	57/72 mm
linearity +-2% peak-peak, Z / XY +- 10% peak-peak, Z / XY	30 / 43 mm 40 / 52 mm
inductance	=< 70 μ H
resistance	=< 1.6 Ω
rise time, 0-40A, 120V	< 150 μ s
cooling	water
maximum current tested	50 A

Table 7.9. B-AFPA-40 Dip Switch Setting

Gradient system	Gradient	C Dip Switch	R Dip Switch	C	R
Mini0.5	X	0001 1100	0000 1000	4.7 nF	15 k Ω
Mini0.5	Y	0001 0000	0000 1000	2.2 nF	15 k Ω
Mini0.5	Z	0001 0000	0001 0000	2.2 nF	20 k Ω


Receiver



Performance of each element in the receiver chain

$$NF_{\text{total}} = NF_{\text{preamp}} + \frac{(NF_{\text{cable}} - 1)}{G_{\text{preamp}}} + \frac{NF_{\text{amp}}}{G_{\text{preamp}} G_{\text{cable}}} + \frac{NF_{\text{attenuator}}}{G_{\text{preamp}} G_{\text{cable}} G_{\text{amp}}} \quad (7.4)$$

Using the numbers in Figure 7.3, the overall NF is 0.51 dB, which is very close to the NF of the first-stage amplifier. Indeed, eqn (7.4) shows that the critical component in terms of the overall NF of the receiver is the first stage,








Figure 7.4 Three variations of T/R switch based on a quarter-wavelength transmission line (a) and its lumped element pi-section equivalent (b) and (c). In (c) two pi-sections are used for additional isolation between transmitter and receiver. Capacitors C_b are very high value to provide a DC-block for the diode-driving voltage, and present essentially a short circuit at RF frequencies.

Resolution of analogue-to-digital converters


3-bit




4-bit



5-bit



Analogue-to-digital converters







Figure 7.7 Block diagram of a FET-based preamplifier connected on the input side to an RF coil impedance matched to 50 Ω, and on the output side to a 50 Ω cable.







Figure 7.7 Block diagram of a FET-based preamplifier connected on the input side to an RF coil impedance matched to 50 Ω , and on the output side to a 50 Ω cable.

Overall control system which integrates everything together




PC (Windows/Linux)


Various sequence descriptions



Red Pitaya SDRlab



Interface with gradient amplifier (optional) and gradient coils



Four channel gradient 16 bit DAC board

Internal power stage which can deliver ± 10 Amps on each channel

Can also output ± 12 volts output for external gradient power amplifier



Clock speed 40 MHz – so if all 4 channels are used in parallel, clocks at 100,000 samples per second, so 10 microseconds per sample


Current is continuously monitored which can be used to calibrate the system.




Compact and mobile scanners are changing workflows in radiology, where MRI comes to the patient rather than the patient to the MRI. In this project we will establish the necessary metrological framework for the harmonised development of clinical low-field MRI.

[Read more](#)





a) Pulseseq Sequence



b) Unrolled Sequence

